


CERTIFICATE OF ANALYSIS

* FOR QUALITY ASSURANCE PURPOSES. NOT A CALIFORNIA COMPLIANCE CERTIFICATE.

THIS IS A REVISED COA ISSUED JUL 08, 2021

SAMPLE: MUSCLE RUB PLUS (+) (TOPICAL) // CLIENT: QUANTA INC. // BATCH: **PASS**



BATCH NO.: AMQP0001
 MATRIX: TOPICAL
 CATEGORY: OTHER
 SAMPLE ID: BCL-210626-006
 COLLECTED ON: JUN 26, 2021
 RECEIVED ON: JUN 26, 2021
 BATCH/SAMPLE SIZE: 28.3 G / 28.3 G
 PACKAGE SIZE: 28.3 G

CANNABINOID OVERVIEW

TOTAL THC:	0 %
TOTAL CBD:	0.54875 %
TOTAL CANNABINOIDS:	1.10653 %

DISTRIBUTOR INFO

DISTRIBUTOR

QUANTA INC.
 632 S GLENWOOD PL
 BURBANK, CA 91506

LICENSE

QUANTA INC.
 ADULT-USE - RETAILER LICENSE



BCL-03: CANNABINOID POTENCY BY HPLC-UV // JUN 29, 2021

TERPENOID TESTING BY HS-GC/FID

ANALYTE	AMT	AMT	LOD/LOQ (mg/g)	PASS/FAIL
CBC	0.00416 %	0.0416 mg/g	0.00102/0.02317	N/A
CBCA	0.00984 %	0.0984 mg/g	0.00292/0.02317	N/A
CBD	0.49291 %	4.9291 mg/g	0.00144/0.02317	N/A
CBDA	0.06367 %	0.6367 mg/g	0.00074/0.02317	N/A
CBDV	0.00297 %	0.0297 mg/g	0.00111/0.02317	N/A
CBDVA	ND	ND	0.00074/0.02317	N/A
CBG	ND	ND	0.00134/0.02317	N/A
CBGA	< LOQ	< LOQ	0.00093/0.02317	N/A
CBN	ND	ND	0.00074/0.02317	N/A
Δ ⁸ -THC	ND	ND	0.00269/0.02317	N/A
Δ ⁹ -THC	ND	ND	0.00195/0.02317	N/A
THCA	ND	ND	0.00139/0.02317	N/A
THCV	0.33285 %	3.3285 mg/g	0.00153/0.02317	N/A
THCVA	0.20013 %	2.0013 mg/g	0.00107/0.02317	N/A
TOT THC **	ND	ND		N/A
TOT CBD **	0.54875 %	5.4875 mg/g		N/A
CBD/PKG	139.49 mg			N/A
Δ ⁹ -THC/PKG	ND			PASS

ANALYTE	AMT	AMT	LOD/LOQ ()	PASS/FAIL
TOTAL TERPENES	ND	ND		N/A
α-PINENE	NT	NT		N/A
CAMPHERE	NT	NT		N/A
β-PINENE	NT	NT		N/A
β-MYRCENE	NT	NT		N/A
Δ ³ -CARENE	NT	NT		N/A
α-TERPINENE	NT	NT		N/A
D-LIMONENE	NT	NT		N/A
P-CYMENE	NT	NT		N/A
α-OCIMENE	NT	NT		N/A
β-OCIMENE	NT	NT		N/A
γ-TERPINENE	NT	NT		N/A
TERPINOLENE	NT	NT		N/A
LINALOOL	NT	NT		N/A
ISOPULEGOL	NT	NT		N/A
GERANIOL	NT	NT		N/A
β-CARYOPHYLLENE	NT	NT		N/A
α-HUMULENE	NT	NT		N/A
CIS-NEROLIDOL	NT	NT		N/A
TRANS-NEROLIDOL	NT	NT		N/A
GUAIOL	NT	NT		N/A
α-BISABOLOL	NT	NT		N/A
1,8-CINEOLE *	NT	NT		N/A

** TOTAL THC = (THCA X 0.877) + THC
 ** TOTAL CBD = (CBDA X 0.877) + CBD

* BEYOND SCOPE OF ACCREDITATION



RESULTS CERTIFIED BY: RON BROOKS, M.S.
 LABORATORY DIRECTOR, BELCOSTA LABS
 JUL 08, 2021

